

INVOTEC® FINAL™

ABSORBABLE NASAL DRESSING & PACK

Invotec International Final™ is formulated to be used as a nasal dressing to prevent adhesions and to stop bleeding.

- Prevents nasal adhesions.
- Performs fast and effective hemostasis.
- Contributes to advanced wound healing.
- Liquefies in 3 days, absorbs completely in 7-9 days.

THE FINAL™ ANSWER TO NASAL DRESSINGS...

sixty years of proven clinical effectiveness and safety, just got better!

INVOTEC® FINAL™

ABSORBABLE NASAL DRESSING & PACK



Formulated for use after nasal surgery to prevent adhesions and to treat epistaxis.

Superior Hemostasis¹

Collagen gelatin induces the adherence and aggregation of platelets, which initiates hemo-coagulation and absorbs greater than 45 times its own weight. Stops bleeding one (1) minute after insertion.

Direct Wound Healing²

Provides a superior wound-healing environment, not just by maintaining a moist environment, but also due to its material composition contributing directly in achieving more rapid healing, healthier and better structured wounds.



Absorption & Liquification

Final liquefies in 3 days and absorbs completely within 7-9 days.

Intended Use: The Invotec® Final™ Nasal Dressing and Pack is intended to prevent PostOP nasal adhesions (synechia) and to stop nasal bleeding.

Material Composition: The Invotec® Final™ is made from 100% purified porcine gelatin, especially formulated to be used in ENT procedures.

¹ M. Suckfull MD, Ludwig-Maximilian University, Munich, Germany, Clinical Observations: Oct 2009. Document on file.

² L.Doorn, MSc, Susan Klymowsky MD, The Effect of Gelatin on Wound Healing, Consultative Report, Qserve Amsterdam, Oct, 2009.

Catalog No.	Description	Qty/Box
GS805004	Bilateral, Size 80mm x 50mm x 4mm	10
GS802504	Single Side (Small) 80mm x 25mm x 4mm	10

